



## Generate publication-quality western blot data with TrueBlot®

現在ご使用のウェスタンブロットにTrueBlot™を使用するのは簡単です。現在お使いのHRP 標識二次抗体をマウス、ウサギ、ヤギ、ヒツジ TrueBlot™ HRP Ig に変えるだけです。TrueBlot™はマウス、ウサギ、ヤギ、ヒツジIgGの native なジスルフィド型のみを選択的に検出します。  
 - 還元された抗体重鎖と抗体軽鎖は検出されません

免疫沈降、ウェスタンブロットに用いることができるTrueBlot™製品:

TrueBlot™ Immunoprecipitation Beads  
 TrueBlot™ HRP IgG  
 TrueBlot™ Sets  
 TrueBlot™ Western Blot Kits

## TrueBlot® Ordering Information

TrueBlot® Immunoprecipitation (IP) Beads	Size	Cat. No.	希望小売価格
TrueBlot® anti-Mouse Ig IP Beads	2.5 ml	00-8811-25	¥39,100
TrueBlot® anti-Rabbit Ig IP Beads	2.5 ml	00-8800-25	¥39,100
TrueBlot® anti-Goat Ig IP Beads	2.5 ml	00-8844-25	¥39,100
Mouse TrueBlot® (HRP) Products	Size	Cat. No.	
Mouse TrueBlot® ULTRA: HRP anti-mouse IgG	20 µl (1000X)	18-8817-30	¥9,200
	50 µl (1000X)	18-8817-31	¥19,600
	200 µl (1000X)	18-8817-33	¥42,600
Mouse TrueBlot® ULTRA: Biotin anti-mouse IgG	25 µg	13-8817-80	¥9,200
	100 µg	13-8817-82	¥34,500
Mouse TrueBlot® Set (includes 50 µl of Mouse TrueBlot® ULTRA & 2.5 ml IP Beads)	50 µl (1000X)	88-7788-31	¥56,400
Mouse TrueBlot® Western Blot Kit (includes 50 µl of Mouse TrueBlot® ULTRA, substrate solutions & blocking buffers)	50 µl (1000X)	88-8887-31	¥54,100
Rabbit TrueBlot® (HRP) Products	Size	Cat. No.	
Rabbit TrueBlot®: HRP anti-rabbit IgG	50 µl (1000X)	18-8816-31	¥21,900
	200 µl (1000X)	18-8816-33	¥44,900
Rabbit TrueBlot® Set (includes 50 µl of Rabbit TrueBlot® & 2.5 ml IP Beads)	50 µl (1000X)	88-1688-31	¥58,700
Rabbit TrueBlot® Western Blot Kit (includes 50 µl of Rabbit TrueBlot®, substrate solutions & blocking buffers)	50 µl (1000X)	88-8886-31	¥56,400
Goat TrueBlot® (HRP) Products	Size	Cat. No.	
Goat TrueBlot®: HRP anti-goat IgG	50 µl (1000X)	18-8814-31	¥21,900
	200 µl (1000X)	18-8814-33	¥44,900
Goat TrueBlot® Set (includes 50 µl of Goat TrueBlot® & 2.5 ml IP Beads)	50 µl (1000X)	88-1488-31	¥58,700
Goat TrueBlot® Western Blot Kit (includes 50 µl of Goat TrueBlot®, substrate solutions & blocking buffers)	50 µl (1000X)	88-8884-31	¥56,400
Sheep TrueBlot® (HRP) Products	Size	Cat. No.	
Sheep TrueBlot®: HRP anti-sheep IgG	50 µl (1000X)	18-8815-31	¥21,900
	200 µl (1000X)	18-8815-33	¥44,900

### ■ 輸入販売元



ベイ バイオサイエンス株式会社  
 テクニカルサポート TEL :078-304-5881  
 Home Page: <http://www.baybio.co.jp>

### ■ 代理店

### ■ 東日本 販売元 ■ 中部

ノベルサイエンス株式会社  
 NOVEL SCIENCE CO., Ltd.

〒113-0033 東京都文京区本郷7丁目2番1号  
 TEL (03)5842-2901 FAX (03)5842-2905

### ■ 西日本 販売元

和研薬株式会社  
 WAKENYAKU CO., LTD.

本社 〒606-8171 京都市左京区一乗寺西水干町17番地  
 西日本営業部 〒525-0029 滋賀県草津市下笠町945番地1  
 TEL (077)568-2911/FAX (077)568-2944

本社営業部 (075)721-8111 京阪奈営業所 (0774)65-2521  
 滋賀営業所 (077)568-2481 神戸営業所 (078)306-6060  
 大阪営業所 (072)749-5200 <http://www.wakenyaku.co.jp>